

**Medtrica Solutions, Ltd**  
**4001 Inglewood Ave; 101-128**  
**Redondo Beach, CA 90278**  
**P: 877-703-1079/F: 310-263-9098**  
[www.medtrica.com](http://www.medtrica.com)



## **ESW905200**

Flat Sponge: w/ 30 ml of pre-mixed Enzymatic Detergent SOLUTION,  
Individually Sealed

### **Directions for Use:**

1. Open package using the tear notches and remove the sponge
2. Use the sponge to wipe down the exterior of the endoscope
3. Discard the sponge using approved methods for materials that have been in contact with possible infectious waste

### **Product Details:**

- For precleaning of instruments after a procedure. Low foaming formula
- Enzymatic Detergent concentrate contains 3 enzymes - Protease, Amylase and Lipase, PLUS Surfactant Detergents to break down proteinaceous bioburden during the preclean of surgical instruments
- The water that is mixed with enzymatic detergent concentrate goes through a multi-stage cleaning process in the following order:  
1) Particulate Filtration; 2) Charcoal Filtration; 3) Semi Permeable Membrane (reverse osmosis); 4) UV Light
- Flat Sponge size: 5"L x 4"W x .5"H (127mm x 102mm x 12.5mm)
- Green Friendly, High Density, 100% Medical Grade foam, No Flaking or Scratching
- Biodegradable, non-toxic, odorless and non-irritating to the eyes, bacterio static
- Notched, easy tear open packaging
- Bulk Packaged - 200 ct /1 case or Dispenser Boxes - 4 boxes of 50 - (200 ct / 1 case)