



Bedside Cleaning Add Water Products





MEDTRICA

D-Zyme™ Concentrate Sponge

Description	Code	
<p>For optimal rinsing and pre-cleaning of instruments at bedside</p> <p>Activate concentrate with up to 500 ml of water</p> <p>Delivers Low Foaming D-Zyme™ enzymatic detergents containing 3 enzymes - Protease, Amylase and Lipase, PLUS Surfactant Detergents to break down proteinaceous bioburden during the pre-clean of surgical instruments</p> <ul style="list-style-type: none"> Green Friendly, 100% Medical Grade foam, No Flaking or Scratching Biodegradable, non-toxic, odorless and non-irritating to the eyes, bacteriostatic <p>Options:</p> <ul style="list-style-type: none"> D-Zyme™ Concentrate Flat - Add to Water - 200 ct / 1 case D-Zyme™ Concentrate Tubular - Add to Water - 100 ct / 1 case 	<p>EWC805200</p> <p>EWC804100</p>	 <p>D-Zyme™ Concentrate Flat</p>  <p>D-Zyme™ Concentrate Tubular</p>

D-Zyme™ Concentrate with 500 ml Basin

<p>Add up to 500ml of water to the clear plastic basin, then soak and squeeze the sponge to activate D-Zyme™ Enzymatic Concentrate</p> <p>Options:</p> <ul style="list-style-type: none"> D-Zyme™ Concentrate Flat w/ 500 ml clear Basin - 100 ct / 1 case D-Zyme™ Concentrate Tubular w/ 500 ml clear Basin - 50 ct / 1 case <p><i>Also available:</i></p> <p>**D-Zyme™ Concentrate Flat in Master Bag & 500ml Basin - Add Water - 50 ct / 1 case</p> <p>**D-Zyme™ Concentrate Tubular in Master Bag & 500ml Basin - Add Water - 24 ct / 1 case</p>	<p>ETF500100</p> <p>ETT500050</p>	 <p>D-Zyme™ Flat Concentrate & 500 ml Basin</p>  <p>D-Zyme™ Tubular Concentrate & 500 ml Basin</p>
---	-----------------------------------	--

Manufactured by:

MEDTRICA SOLUTIONS LTD.
104A-6350 204th Street
Langley, BC Canada V2Y 2V1
t: 604-618-2291



info@medtrica.com
www.medtrica.com

MEDTRICA SOLUTIONS LTD

4001 Inglewood Ave
101-128
Redondo Beach, CA 90278
t: 877-703-1079