

Bedside Cleaning Ready to Use



MEDTRICA

D-ZYME™ SPONGES - PRE-SOAKED w/ 30ml of D-ZYME™ SOLUTION

Description	Code	
<p>Flat D-Zyme™ sponge: for pre-cleaning instruments and Surfaces/ Multi-Use</p> <p>Tubular D-Zyme™ sponge: for pre-cleaning Endoscopes - and as a bumper guard during transport</p> <p>D-Zyme™ sponges absorb more cleaning solution than other sponges/ paper based wipes - Allows for the safe cleaning of sensitive instruments and fixtures, as well as general surfaces</p> <p>Delivers 30 ml of Low Foaming D-Zyme™ enzymatic detergents solution containing 3 enzymes - Protease, Amylase and Lipase, PLUS Surfactant Detergents to break down proteinaceous bioburden during the pre-clean of surgical instruments</p> <ul style="list-style-type: none"> • Green Friendly, 100% Medical Grade foam, No Flaking or Scratching • Biodegradable, non-toxic, odorless and non-irritating to the eyes, bacteriostatic <p>The water that is mixed with D-Zyme™ enzymatic detergents concentrate to make the pre-mixed solution goes through a mult-stage cleaning process in the following order: 1) Particulate Filtration; 2) Charcoal Filtration; 3) Semi Permeable Membrane (reverse osmosis); 4) UV Light</p>		 <p>D-Zyme™ Flat</p>
<p>Options:</p> <p>Flat D-Zyme™ - 30 ml</p> <ul style="list-style-type: none"> • Bulk - 200 ct / 1 case • Dispenser Box - 200 ct / 1 case - (4 boxes of 50) <p>Tubular D-Zyme™ - 30 ml</p> <ul style="list-style-type: none"> • Bulk - 100 ct / 1 case • Dispenser Box - 100 ct / 1 case - (4 boxes of 25) 	<p>ESW905200 ESW905204</p> <p>ESW904100 ESW904104</p>	 <p>D-Zyme™ Tubular</p>  <p>Dispenser Box</p>

Manufactured by:

MEDTRICA SOLUTIONS LTD.
104A-6350 204th Street
Langley, BC Canada V2Y 2V1
t: 604-618-2291



info@medtrica.com
www.medtrica.com

MEDTRICA SOLUTIONS LTD

4001 Inglewood Ave
101-128
Redondo Beach, CA 90278
t: 877-703-1079