

OR / Surgery Table Pads



MEDTRICA



OR Table Pads Features:

CHOICE OF OR TABLE PADS CONSTRUCTION:

- Single Layer High Resiliency Foam
- Dual Layer - HD Base Foam With a 5.0 lb Memory Foam Top Layer

COVERS:

All covers are 100% medical grade fabric - antifungal, antibacterial, latex free and completely fluid-proof. Choices are a durable Vinyl or a Pressure Reducing stretchable Dartex™ fabric. All covers are fluid proof, anti-fungal, anti-bacterial and Latex free.

FOAM CORES:

100% medical grade Green Friendly non-recycled, non-allergenic, HD polyurethane and Viscoelastic (memory) foam without fillers or additives.

CLEANING:

OR Table Pad is disinfected and wiped clean using standard cleaning products.

GUARANTEE:

A full non prorated guarantee on both the cover and the core for all OR Table Pads

Adhesive Velcro included to keep the mattress secure to the OR Table top.

Custom Sizing and Construction Available at No Extra Cost.

**MORE DURABLE,
MORE VERSATILE**

Medtrica OR Table Pads provide greater comfort and durability because of unsurpassed quality. Manufactured in the USA and Canada - built to last.

100% medical grade Green Friendly foam of varying Densities and IFD are used for the construction of each foam layer.

Dartex™ Stretchable Pressure Reducing Fabrics are the world leader in Healthcare cover fabrics.

Dartex™ fabrics are engineered to provide the best in MVT, pressure management, abrasion and tear resistance while maintaining optimal performance throughout the entire lifespan of the stretcher mattress.

FLAMMABILITY STANDARDS

All materials meet Cal TB 117 and CFR 1632
Option to meet Cal TB 129 & Boston FC IX-11

AVAILABLE SIZES:

All Standard and Custom sizes available to Fit All Makes and Models of OR/Surgery Tables.

All OR Table Pads Meet and Exceed OEM specifications

Individual Pad Sections Available

Armboard Pads Available

Table Straps and Positioners available

A better quality Surgery Table Pad - Safer and more Durable - for a better price

Manufactured by:

MEDTRICA SOLUTIONS LTD.
6290, 204th Street, #104A
Langley, BC. Canada V2Y 2V1
t: 604-618-2291
f: 604-952-4314



info@medtrica.com
www.medtrica.com

MEDTRICA SOLUTIONS LTD

4001 Inglewood Ave
101-128
Redondo Beach, CA 90278
t: 877-703-1079
f: 310-263-9098