


# SolaCare Psychiatric Standard



MEDTRICA



## NON-ALLERGENIC

The entire foam core is constructed of 100% medical grade Green Friendly non-allergenic foam.

## ANTI-MICROBIAL

All components are anti-microbial and assist in prohibiting growth of bacteria and fungi.

## FLAMMABILITY

Meets Federal Standard 16 CFR 1632, Cal. TB 117. Optional barrier to meet Cal. TB 129 and Boston Fire IX-11.

## COVER DESIGN

Fully enclosed cover, No Handles, No Flaps, No Zippers.

## PSYCHIATRIC STANDARD MATTRESS COVER

Durable, snug fitting low denier 70D Nylon is soft, moisture vapor permeable, and resistant to staining with low slip vinyl bottom cover. Will not allow fluids to pass through and can be easily cleaned with standard hospital cleaning products. 100% 70D Nylon, no butyl backing.

## OPTIONAL MATTRESS COVER - DARTEX™

Dartex™ 4 way Stretch Pressure Reducing fabrics are the world leader in healthcare fabrics. Dartex™ fabrics are engineered to provide the best in MVT, pressure management, abrasion and tear resistance while maintaining optimal performance throughout the entire lifespan of the product.

## GUARANTEE

A full non prorated guarantee on all SolaCare Series mattresses instead of the industry standard limited warranty.

- 3 year cover guarantee
- 3 year foam core guarantee

## FUNCTIONAL & DURABLE

*The SolaCare Psychiatric Standard Mattress is specifically engineered for Durability, Support and Patient Safety using the highest quality medical grade materials.*

*Reversible design (& Flippable w/ optional cover) has a Fully Enclosed cover with No Handles, No Flaps, No Zippers. The Psychiatric Standard accommodates patients up to 300 lbs and is Latex Free.*

info@medtrica.com  
www.medtrica.com



**A better quality mattress, more durable, safe and comfortable - for a better price.**