


Safe Room Mattress



MEDTRICA



NON-ALLERGENIC

The entire foam core is constructed of 100% medical grade Green Friendly non-allergenic foam.

ANTI-MICROBIAL

All components are anti-microbial and assist in prohibiting growth of bacteria and fungi.

FLAMMABILITY

Meets Federal Standard 16 CFR 1632, Cal. TB 117. Optional barrier to meet Cal. TB 129 and Boston Fire IX-11.

COVER DESIGN

Fully enclosed, No Handles, No Flaps, No Zippers. Tear resistant w/ Rip-stop weave.

MATTRESS COVER

Durable, snug fitting Heavy Duty, Tear resistant Vinyl cover w/ Rip-stop weave is moisture vapor permeable, and resistant to staining. Will not allow fluids to pass through and can be easily cleaned with standard hospital cleaning products.

GUARANTEE

A full non prorated guarantee on all SolaCare Series mattresses instead of the industry standard limited warranty.

- 2 year cover guarantee
- 2 year foam core guarantee

FUNCTIONAL & DURABLE

The Safe Room Mattress is specifically engineered for Durability, Support and Patient Safety using the highest quality medical grade materials.

Reversible & Flippable design has a Fully Enclosed cover with No Handles, No Flaps, No Zippers. Heavy Duty Vinyl cover is Tear resistant w/ Rip-stop weave.

The Safe Room Mattress accommodates up to 450 lbs and is Latex Free.

info@medtrica.com
www.medtrica.com



A better quality mattress, more durable, safe and comfortable - for a better price.