

# SolaCare Memory Flex 3.0 XT



MEDTRICA

## Multi Zoned & Layered Foam Support

### Anti-Shear Liner option:

Dissipates Heat and Reduces Friction

### Pressure Reducing Top Foam Support

3.0 lb density Viscoelastic (memory foam) covering the entire patient surface of the mattress.

### Perimeter Foam Support

Constructed of 2.0lb density and IFD 86-94.

### Head section

2.2lb density and IFD 23-31.

### Seat section

Ultracel™ 2.5lb density and IFD 24-32.

### Fluid Proof Cover Options:

- Dartex™ 4 way stretch
- 70D Nylon

### Base Foam Support

Constructed of 2.2lb density and IFD 51-59. Provides a firm support base that prevents "bottoming out" even as the head section is fully elevated.

### Calf section

2.2lb density and IFD 23-31.

### Heel section

1.6lb density and IFD 16-24.

### AVAILABLE SIZES:

- 80" x 54" x 6"
- 80" x 48" x 6"
- 80" x 42" x 6"
- 80" x 39" x 6"
- 75" x 54" x 6"
- 75" x 39" x 6"

### NON-ALLERGENIC

The entire foam core is constructed of 100% medical grade Green Friendly non-allergenic foam.

### ANTI-MICROBIAL

All components are anti-microbial and assist in prohibiting growth of bacteria and fungi.

### FLAMMABILITY

Meets Federal Standard 16 CFR 1632, Cal. TB 117. Optional barrier to meet Cal. TB 129 and Boston Fire IX-11.

### CLAM-SHELL DESIGN COVER

Clam-shell design cover w/Double Folded Seams and 4 Reinforced Handles.

### STANDARD MATTRESS COVER

Durable, snug fitting low denier 70D Nylon is soft, moisture vapor permeable, and resistant to staining with low slip vinyl bottom cover. Will not allow fluids to pass through and can be easily cleaned with standard hospital cleaning products. 100% 70D Nylon, no butyl backing.

### OPTIONAL MATTRESS COVER

Dartex™ 4 way Stretchable Pressure Reducing Fabrics are the world leader in Healthcare fabrics. Dartex™ fabrics are engineered to provide the best in MVT, pressure management, abrasion and tear resistance while maintaining optimal performance throughout the entire lifespan of the product.

### GUARANTEE

A full non prorated guarantee on all SolaCare Series mattresses instead of the industry standard limited warranty.

- 5 year cover guarantee
- 5 year foam core guarantee

### SAFETY, COMFORT & DURABLE

The SolaCare Memory Flex 3.0 XT Bariatric Mattress is specifically engineered for lasting optimum performance by using better pressure reducing materials. Comfort, safety and durability is enhanced by the: 3 Layer, 4 Zone design (including the top 3lb memory foam layer), Firm internal foam side rails, Anti-Shear Liner and Dartex™ fabric cover. Angled Heel section and other options available. The Memory Flex 3.0 XT Bariatric accommodates patients up to 1000 lbs and is Latex Free.

info@medtrica.com  
www.medtrica.com



**A better quality mattress, more durable, safe and comfortable - for a better price.**