

Straps/Belts



STRAPS

Straps are made with either a soft durable Nylon, a tear resistant Vinyl or Rubber that are easy to wipe clean. All materials exceed flammability standards Cal TB 117 and CFR 1632. Straps are anti-bacterial, anti-fungal, stain and tear resistant, anti-static and latex free. All straps are radiolucent, non-magnetic and artifact-free.


Strap Description		
<p>O.R. Table Strap w/ Vinyl covered padded Abdominal Panel High strength nylon webbing that is soft and easy to clean. Fully adjustable with double pinch quick release buckles. Size: 96"L x 2"W (5"W abdominal panel), ea Size: 96"L x 2"W (5"W abdominal panel), Package of 4</p>	<p>ORA96000500 ORA96000504</p>	
<p>Rubber OR Table Strap - Amsco style Size: 96"L x 2"W (5"W ab panel), w/ Metal Buckles - ea Size: 132"L x 2"W (5"W ab panel), w/ Metal Buckles -ea Size: 96"L x 2"W w/ Metal buckles or Hook - ea Size: 132"L x 2"W w/ Metal buckles or Hook - ea Size: 96"L x 2"W (5"W ab panel), w/ HD Plastic buckles - ea Size: 132"L x 2"W (5"W ab panel), w/ HD Plastic buckles - ea Size: 96"L x 2"W w/ HD Plastic cam action buckles - ea Size: 132"L x 2"W w/ HD Plastic cam action buckles - ea</p>	<p>ORR96050200 ORR13250200 ORR96000200 ORR13200200 ORP96050200 ORP13250200 ORP96000200 ORP13200200</p>	
<p>Disposable padded abdominal panel wrap for all rubber Amsco style OR Table straps - package of 10</p>	<p>ABR140500</p>	
<p>O.R. Positioning Straps - Nylon Durable soft Nylon strap with hook and loop on the same side. Straps are secured through the rails of the table. Size: 60"L x 2"W, ea Size: 60"L x 2"W, Package of 4 Size: 76"L x 2"W, ea Size: 76"L x 2"W, Package of 4 Size: 132"L x 2"W, ea</p>	<p>ORN60000200 ORN60000204 ORN76000200 ORN76000204 ORN13200200</p>	
<p>O.R. Positioning Straps - Vinyl Durable Vinyl strap with hook and loop on the same side. Straps are secured through the rails of the table. Size: 60"L x 2"W, ea Size: 60"L x 2"W, Package of 4 Size: 76"L x 2"W, ea Size: 76"L x 2"W, Package of 4 Size: 60"L x 4"W, ea Size: 60"L x 4"W, Package of 4 Size: 76"L x 4"W, ea Size: 76"L x 4"W, Package of 4 Size: 60"L x 12"W, ea Size: 60"L x 12"W, Package of 4 Size: 76"L x 12"W, ea Size: 76"L x 12"W, Package of 4 Size: 60"L x 24"W, ea Size: 60"L x 24"W, Package of 4 Size: 76"L x 24"W, ea Size: 76"L x 24"W, Package of 4 Size: 132"L x 2"W, ea</p>	<p>ORV60000200 ORV60000204 ORV76000200 ORV76000204 ORV60000400 ORV60000404 ORV76000400 ORV76000404 ORV76000404 ORV60001200 ORV60001204 ORV60001204 ORV60002400 ORV60002404 ORV76002400 ORV76230404 ORV13200200</p>	

** Custom Sizing and Quantity Discounts Available

*** Full 1 year warranty***

STRAPS

Straps are made with either a soft durable Nylon, a tear resistant Vinyl or Rubber that are easy to wipe clean. All materials exceed flammability standards Cal TB 117 and CFR 1632. Straps are anti-bacterial, anti-fungal, stain and tear resistant, anti-static and latex free. All straps are radiolucent, non-magnetic and artifact-free.

Strap Description		
<p>Medical Imaging Positioning Straps - Vinyl Durable Vinyl strap with hook and loop on the opposite side. Straps are secured through the rails of the table.</p> <p>Size: 96"L x 2"W (5"W abdominal panel), ea Size: 60"L x 2"W, ea Size: 60"L x 2"W, Package of 4 Size: 76"L x 2"W, ea Size: 76"L x 2"W, Package of 4 Size: 60"L x 4"W, ea Size: 60"L x 4"W, Package of 4 Size: 76"L x 4"W, ea Size: 76"L x 4"W, Package of 4 Size: 60"L x 12"W, ea Size: 60"L x 12"W, Package of 4 Size: 76"L x 12"W, ea Size: 76"L x 12"W, Package of 4 Size: 60"L x 24"W, ea Size: 60"L x 24"W, Package of 4 Size: 76"L x 24"W, ea Size: 76"L x 24"W, Package of 4</p>	<p>ORA96000500 MIV60000200 MIV60000204 MIV76000200 MIV76000204 MIV60000400 MIV60000404 MIV76000400 MIV76000404 MIV60001200 MIV60001204 MIV76001200 MIV76001204 MIV60002400 MIV60002404 MIV76002400 MIV76002404</p>	
<p>Arm Position Straps - Nylon Durable soft Nylon strap with hook and loop to easily secure on to the wrists or ankles.</p> <p>Size: 60"L x 2"W (12"L - 'T' section), ea Size: 60"L x 2"W (12"L - 'T' section), Package of 4 Size: 76"L x 2"W (12"L - 'T' section), ea Size: 76"L x 2"W (12"L - 'T' section), Package of 4</p>	<p>TSN60000200 TSN60000204 TSN76000200 TSN76000204</p>	
<p>Arm Board Straps - Vinyl Durable Vinyl strap w/ HD plastic cam action buckle</p> <p>Size: 1"W x 32"L - Package of 10</p>	<p>OAV3300100</p>	
<p>Head Positioning Straps - Nylon Durable soft Nylon strap with hook and loop on the same side for easy positioning. The strap is secured through the side rails. Head straps are non-magnetic and radiolucent.</p> <p>Size: 30"L x 1"W, ea Size: 30"L x 1"W, Package of 4 Size: 30"L x 2"W, ea Size: 30"L x 2"W, Package of 4 Size: 36"L x 1"W, ea Size: 36"L x 1"W, Package of 4 Size: 36"L x 2"W, ea Size: 36"L x 2"W, Package of 4</p>	<p>HSN30000100 HSN30000104 HSN30000200 HSN30000204 HSN36000100 HSN36000104 HSN36000200 HSN36000204</p>	

** Custom Sizing and Quantity Discounts Available

* Full 1 year warranty*

Manufactured by:

MEDTRICA SOLUTIONS LTD.
6350 204th Street, #104A
Langley, BC. Canada V2Y 2V1
t: 604-618-2291
f: 604-952-4314



info@medtrica.com
www.medtrica.com

MEDTRICA SOLUTIONS LTD

4001 Inglewood Ave
101-128
Redondo Beach, CA 90278
t: 877-703-1079
f: 310-263-9098

STRAPS

An all-in one Velcro™ that is hook on one side and loop on the opposite side. One-Wrap® is a both an economical solution to any application that requires positioning or securing. Can be cut into smaller lengths

*** Full 1 year warranty***

One-Wrap® Positioning Strap	Product Code	
Size: 60"L x 0.5"W, ea	OWV60007500	
Size: 60"L x 0.5"W, Package of 4	OWV60007504	
Size: 76"L x 0.5"W, ea	OWV76007500	
Size: 76"L x 0.5"W, Package of 4	OWV76007504	
Size: 60"L x 1"W, ea	OWV60000100	
Size: 60"L x 1"W, Package of 4	OWV60000104	
Size: 76"L x 1"W, ea	OWV76000100	
Size: 76"L x 1"W, Package of 4	OWV76000104	
Size: 60"L x 2"W, ea	OWV60000200	
Size: 60"L x 2"W, Package of 4	OWV60000204	
Size: 76"L x 2"W, ea	OWV76000200	
Size: 76"L x 2"W, Package of 4	OWV76000204	
Size: 60"L x 4"W, ea	OWV60000400	
Size: 60"L x 4"W, Package of 4	OWV60000404	
Size: 76"L x 4"W, ea	OWV76000400	
Size: 76"L x 4"W, Package of 4	OWV76000404	
** Custom Sizing and Quantity Discounts Available		

SANDBAGS

All Sandbags are made with double-wall construction. The outer wall consists of a soft durable Vinyl that is tear resistant and easy to wipe clean. All materials exceed flammability standards Cal TB 117 and CFR 1632. Sandbags are anti-bacterial, anti-fungal, stain resistant, anti-static and latex free. All straps are radiolucent, non-magnetic and artifact-free.

Sandbag Description	Product Code	
O.R. / Medical Imaging Sandbags Size: 9"L x 5"W x 2.5"H Size: 11"L x 5"W x 2.5"H Size: 16"L x 6"W x 2.5"H Size: 12"L x 12"W x 2.5"H	SBV090524 SBV110525 SBV160625 SBV121225	
O.R. / Medical Imaging Sandbags w/ Handle Size: 9"L x 5"W x 2.5"H Size: 11"L x 5"W x 2.5"H Size: 16"L x 6"W x 2.5"H Size: 12"L x 12"W x 2.5"H	SHV090525 SHV110525 SHV160625 SHV121225	
Donut Sandbag Size: 9"Dia x 2.5"H	SDV090025	
** Custom Sizing and Quantity Discounts Available		

Manufactured by:

MEDTRICA SOLUTIONS LTD.
 6350 204th Street, #104A
 Langley, BC, Canada V2Y 2V1
 t: 604-618-2291
 f: 604-952-4314



info@medtrica.com
 www.medtrica.com

MEDTRICA SOLUTIONS LTD

4001 Inglewood Ave
 101-128
 Redondo Beach, CA 90278
 t: 877-703-1079
 f: 310-263-9098