


# UltraFlex Stretcher Mattress



MEDTRICA



## GUARANTEE:

A full non-prorated guarantee on both the cover and the core for all stretcher mattresses.

Adhesive Velcro included to keep the mattress secure to the stretcher.

## FLAMMABILITY

All materials meet Cal TB117 and CFR1632. Option to meet Cal TB129 & Boston FC IX 11.

All Stretcher mattresses meet and exceed OEM specifications.

## CHOICE OF MATTRESS CONSTRUCTION:

- 4" or 5" Thick; All standard and custom Lengths & Widths
- Tapered or Square corners.

## COVERS:

All covers are 100% medical grade fabric-anti-fungal, anti-bacterial, latex free and completely fluid-proof. Choices are a durable Vinyl or a Pressure Reducing stretchable Dartex™ fabric. UltraFlex with Pressure Reducing stretchable Dartex™ fabric includes Anti-Shear liner.

## DARTEX™ 4 WAY STRETCHABLE PRESSURE REDUCING FABRICS:

Dartex™ 4 way Stretchable Pressure Reducing Fabrics are the world leader in Healthcare cover fabrics. Dartex™ fabrics are engineered to provide the best in MVT, pressure management, abrasion and tear resistance while maintaining optimal performance throughout the entire lifespan of the stretcher mattress.

## FOAM CORES:

100% medical grade Green Friendly non-recycled, non-allergenic, HD polyurethane and Viscoelastic (memory) foam without fillers or additives.

## CLEANING:

Stretcher mattresses are disinfected and wiped clean using standard cleaning products.

## MORE DURABLE & VERSATILE

Medtrica stretcher mattresses provide greater comfort & durability because of unsurpassed Quality.

100% medical grade Green Friendly foam of varying Densities & IFD are used for the construction of each foam layer.

Manufactured in the USA & Canada - Engineered to Last.

info@medtrica.com  
www.medtrica.com



*A better quality stretcher mattress, more durable - for a better price.*